

DispensMate

Bottle-Top Dispenser

Easy-to-use volume adjustment

Smooth piston design provides effortless dispensing.

Closure cap protects against reagent contact easy to open/close without taking gloves off.



Adapters allow fitting to most common bottles.

Filling tube adjusts easily for different size bottles.

Volume Range (ml)	Graduation (ml)	AST		CVS	
		%	µl	%	µl
0.5-5	0.1	0.5	25	0.1	5
1.0-10	0.2	0.5	50	0.1	10
2.5-25	0.5	0.5	125	0.1	25
5.0-50	1	0.5	250	0.1	50

Description:

- Excellent chemical resistance, components are made of PTFE, FEP, BSG, PP.
- Fully autoclavable at 121°C.
- Four ranges of bottle-top dispenser covering a volume range from 0.5mL to 50mL.
- Easy for cleaning and maintenance.
- The optional flexible discharge tube with safety handle permits fast and precise dispensing.
- Vapor pressure max. 500mbar, viscosity max. 500mm²/s.
- Temperature max. 40°C, density max. 2.2g/cm³.
- DispensMate is supplied with S40, GL32, GL38, GL25, GL28.

Accessories:

GL45, Adapter, 45mm	S40, Adapter, 45/40mm	GL32, Adapter, 45/32mm
GL38, Adapter, 45/38mm	GL25, Adapter, 32/25mm	GL28, Adapter, 32/28mm
Reagent Bottle (Brown, 1L. Optional)	Filling tube (m)	Discharge tube (m)

