

Huanghua Faithful Instrument Co., Ltd.

www.faithful.cc

F-1B Series Precision Electronic Balance



Characteristic:

The LED display are suit for most use place.F-B series is widely used in industrial applications,laboratory and university because of it made of stable quality.We can use the under weighing functions to test the density.

Main Technical Parameter:

Model	F6001B	F10001B	F20001B	F30001B	F50001BF
Capacity(g)	610	1100	2100	3100	5100
Resolution(0.1				
g)					
Min	0.4				
weighing(g)					
Stable time	≤2S				
Pan size	Ø130mm			160*160mm	
Operation	5-25℃				
temp					
Repeat	0.2g				
ability					
Liner	0.3g				
Cal. weight	200g	500g	1000g	2000g	
Option	RS232C & RS485 \ Printer				
Packing size	Inside:310*230*130mm Outside:660*470*320mm				
G.weight	2.0/22kg				