

Pharmacy Refrigerator

FSF-5V1006

Pharmacy refrigerators are specially designed to store medicines, vaccines, reagents and biomedical products.



Features

- Adjustable temperature range 2-8°C
- Microprocessor control, digital display, and temperature adjustment with an increment of 0.1°C
Large digital display screen
- Temperature variation within $\pm 3^{\circ}\text{C}$
- Permanently lubricated cooling fan for safety and longevity
- High efficiency specialty compressor with known field reliability
- Forced- air cooling with optimized air distribution system designed

Safety

- Build-in back-up battery to display temperature and to operate audible and visual alarm system for up to 8 hours without AC power.
- Visual and audible alarm system: High/low temperature, thermostat failure alarm, power failure alarm, door ajar.
- Power failure protection: turn on delay of the cooling system when power failure, restart delay protection.
- Keyboard lock, password protection to avoid set parameter randomly.
- Wide available voltage range: 187V~242V.

Humanized Design

- 4 casters with stopper for easy moving and fixation.
- 10 adjustable shelves for easy storing.
- Safety lock to prevent unauthorized access .
- Auto on/off LED light with manual switch on control panel.

Reliable control system Inside temperature 2~8°C

Forced-air cooling system Frost-free

Service-free drainage system

Wide voltage range 187V~242V available

5 visual and audible alarms

Performance	
Temperature Range	+2~+8°C
Ambient Temperature	10~32°C
Controller	Microprocessor
Display	Digital Display

Electrical Data	
Power Supply(V/Hz)	220/50~60; 115/60
Rated Power(W)	425
Electrical Current(A)	2.4
Power Consumption	6.23 kWh/24h

Dimensions	
Capacity	1006L
Interior Size(W*D*H)	1100*684*1325mm
Exterior Size(W*D*H)	1240*820*1900mm
Package Size(W*D*H)	1335*945*2000mm
N.t./G.t.(approx)	189/213kg
LoadingQty(20'/40'/40H)	8/18/18

Refrigeration System	
Compressor	1
Refrigerant	R600a
Cooling Type	Forced-air cooling
Defrost Type	Frost free

Alarm	
High/Low Temperature	Y
Power Failure	Y
Door ajar	Y
Thermostat failure	Y

Accessories	
Caster	4
Lamp	LED
Door	2/Heating glass door with automatic rebound
Shelf	10
Test hole	1/25mm
Temperature recorder	Optional