

# Orbital Shaker

## FSK-O330-M Product Description:

- Max. load capacity 3kg;
- Independent LCD display for speed and time;
- Speed and time value can be altered during running;
- Maintenance-free brushless DC motor;
- Foot screws height can be adjusted to make sure stable running;
- Continuous or timing mode optional.



Orbital Diameter	10mm
Max. Load capacity	3kg
Motor type	Brushless DC Motor
Motor Input Power	20W
Motor Output Power	15W
Speed Range	70~400rpm
Speed Display	LCD
Timer Display	LCD
Time Setting Range	1~1199min
Voltage	100-240V, 50/60Hz
Power	25W
Weight	9.8kg
Permissible Ambient temperature and humidity	5~40°C, 80%RH
Dimension(WxDxH)	420×370×100mm
Protection class	IP21

