


Pipette Filler

Levo ME pipette filler has an excellent feature and ergonomic design, suitable for almost all commonly used glass or plastic pipettes.

Levo Me Product Description:

- Enables single-handed operation with minimum effort.
- Compatible with most of the plastic and glass pipettes from 0.1-100mL.
- Easily adjustable speed control.
- 0.45 μ m replaceable hydrophobic filter.
- Stable and reliable performance, elegant appearance.
- High capacity Li-ion battery enable long operation time.



Specifications	Levo Me
Aspiration Speed	25ml<7s
Dispensing Speed	Motor/Gravity
Battery	Lithium-ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2-3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45 μ m Hydrophobic
Weight	200g