

Stainless Steel Water Distiller



Uses:

Stainless steel electric distiller uses electric heating of drinking tap water to produce pure water by distillation, which can be used in medical and sanitary, pharmaceutical, electroplating and laboratory departments. This product is easy to use, simple to operate, with more characteristics of high quality imported stainless steel materials, through special treatment, carefully processed. This not only fully ensures the quality of distilled water, but also greatly improves the service life.

Characteristics:

1. Overall stainless steel design
2. Indicator light shows working status
- 3, convenient operation, reliable performance

Main Parameter:

Model	HS.Z11.5L	HS.Z11.10L	HS.Z11.20L
Output	5L/H	10L/H	20L/H
Material	Stainless steel	Stainless steel	Stainless steel
Input	220V	380V	380V
Heating Power	4.5KW	7.5KW	15KW
NW.	4KG	5KG	8KG
Packing Size(mm)	390×280×750	390×280×800	500×400×970