

JA-N Series Internal Calibration Analytical Balance**Characteristics:**

The visual appearance of the product is blended into the simple concept, stylish and generous. Specially makes the fine pouring cast-aluminium outer covering, so as to enhance the balance of the anti-static anti-jamming capacity. The panel with five operational buttons, sliding glass windscreen provides big space, insightful vision and operates easily. We use high precision magnet sensor to keep the weighing more fast and stable.

Main Technical Parameter:

| MODEL | JA103N | JA203N | JA303N | JA403N | JA503N |
|--------------------------|--|--------|--------|--------|--------|
| Capacity (g) | 100g | 200g | 300g | 400g | 500g |
| Price(USD) | 100g | 200g | 200g | 200g | 200g |
| Resolution(g) | 0.001g | | | | |
| Operating Temp | 5°C-35°C | | | | |
| Repeat ability/Std. Dev. | ±3d | | | | |
| Linearity | ±3d | | | | |
| Weighing Time | ≤ 4s | | | | |
| Sensitivity Drift | ±3d | | | | |
| Corner | ±3d | | | | |
| Weighing Mode | g, oz, ct... | | | | |
| Pan Size | Φ80mm | | | | |
| Key Name | OFF, PRINT, →0←, CAL, ON | | | | |
| Host Dimension | 365X205X350mm | | | | |
| Operation Dimension | 165X170X230mm | | | | |
| Packing Size | Inside:415×285×430mm Outside:475×315×485mm | | | | |
| Display | LCD | | | | |
| Interface | RS232/RJ45/USB... | | | | |
| Internal Calibration | Automatic time trigger calibration | | | | |
| Function expansion | Under weighing / Remote Display... | | | | |